

## Recombinant Enterokinaase

| Catalog Number | Size       |
|----------------|------------|
| AG132-100U     | 100 Units  |
| AG132-1000U    | 1000 Units |

### *Specifications and Use*

#### **Description**

Recombinant Enterokinase (rEK) is the catalytic subunit of bovine enterokinase, which is expressed by the yeast *Pichia pastoris* and purified to yield a high enzyme activity preparation. rEK recognizes the sequence Asp-Asp-Asp-Asp-Lys and cleaves the peptide bond after the lysine residue. The enzyme can be used to cleave any fusion protein that carries this sequence.

#### **Source**

*Yeast*

#### **Unit Definition**

One unit of rEK is the amount of enzyme that will cleave 20µg of thioredoxin-chloramphenicci acetyl transferase fusion protein containing an enterokinase cleavage site (Asp-Asp-Asp-Asp-Lys) to 90% completion at 37°C in 16 hours under the assay conditions listed below.

#### **Assay Conditions**

Recombinant EK in 50mM Tris-HCl, pH 8.0, 1mM CaCl<sub>2</sub>, 0.1% Tween-20, 20µg of fusion protein, and 1 unit rEK in a 30µl reaction volume incubated at 37°C.

#### **Non-Specific Assay**

A non-specific protease activity assay of rEK was performed using azocasein as substrate. The results show that rEK contains less than background levels of non-specific protease.

#### **Storage Conditions**

rEK in 50mM PBNa, pH 8.0, 0.5M NaCl and 50% glycerol should be stored at -20°C. Guaranteed stable for 3 years when stored properly.

### *Recombinant Enterokinaase*

Recombinant Enterokinase (rEK) is a highly specific serine protease that recognizes the amino acid sequence Asp-Asp-Asp-Asp-Lys and cleaves the peptide bond after the lysine residue.

### **FOR RESEARCH USE ONLY**